

Tablet Disintegration Tester

3 Station Park Out



EDI-3X

Features:

- Complies with USP, IP, Ph. Eur. Specifications
- Bathless heating (saves water, time, electricity)
- Bathless disintegration tester is maintenance-free
- Specially designed efficient heater plates for vessel, heats media quickly
- Basket park-out from beaker
- Unique camless drive
- Print report of test results and validation
- Individual motor and timer with external probe for each basket
- Ergonomically designed and magnetically coupled basket for easy loading
- Serial port for PC connectivity

Specifications:

Model	EDI-3X
Beaker	1000 ml glass beakers (3 nos.)
Baskets	3 nos. each carrying 6 test
Guided disc	Perforated cylindrical disc 6
Display	20 x 4 LCD
Heater	M.S.- 400 watts, S.S. 304 watts
Illumination	Green LED
Power supply	230 V AC, 50 Hz, 490 VA
No. of strokes	30 ±0.5 strokes/min
Stroke height	55 mm ±2 mm
Readable temperature range	10° C to 45° C
Controlled temperature range	15° C to 40° C
Temperature resolution	±0.2° C
Timer setting	74 hrs 59 min. 59 sec.

Ordering information:

Products		
Product code	Description	Quantity
0241600	EDI-3X (Body made of S.S.) -230V	1
0221600	EDI-3X (Body made of M.S.) -230V	1
0211600	EDI-3X (Body made of S.S.) -110V	1
0231600	EDI-3X (Body made of M.S.) -110V	1

Standard accessories		
Part number	Description	Quantity
0214A00001	Six tube basket assembly	3
0216A00001	Guided disc	18
0213A00001	Glass beaker (1000 ml)	3
0213B00002	External temperature probe	3

Optional accessories		
Part number	Description	Quantity
0201A00010	Three tube bolus basket assembly with guided disc	3
0201A00007	Magnetic guided disc	18
0201A00011	Single tube bolus basket assembly with guided disc	3
0201B00053	PC keyboard	1
0201A00020	Printer Epson LX-310 Plus	1
0104A00442	Deconex 12 basic 0.7 kgs (approx. 1000ml bottles)	10