

# Hermann

Laboratory Products



*Your professional manufacturing supplier  
of Laboratory Products*

- MAGNETIC HOTPLATE STIRRERS
- OVERHEAD STIRRERS
- VORTEX MIXERS
- ROTATORS
- TUBE ROLLER
- SHAKERS
- HOMOGENIZERS

### » Tube Roller |

Hermann tube rollers provide a gentle and highly efficient rolling motions. Ideal for mixing blood samples, viscous substances and liquid-solid suspensions.

Los rodillos de tubo de Hermann proporcionan movimientos de balanceo suaves y altamente eficientes. Ideal para mezclar muestras de sangre, sustancias viscosas y suspensiones líquido-sólidas.

## HTR-NZ-R12-A

### Tube Roller



#### Performance Highlights

- Speed range of 0~70rpm
- Six roller design for gentle rolling and easy to clean

#### Características Destacadas

- Rango de velocidad de 0 ~ 70 rpm
- Diseño de seis rodillos para rodar suavemente, y fácil de limpiar

Specifications	HTR-NZ-R12-A
Amplitude	24 mm
Voltage	100~240V,50/60Hz
Motor type	DC motor
Timer display	-
Power	25W
Timer	-
Speed range	0~70 rpm
Speed display	Scale
Operation Mode	Continuous
Dimension (DxWxH)	460x260x120mm
Permissible ambient temperature and humidity	5~40 °C, 80%
Protection class	IP21
Rocking Motion	Rocking and Rolling
Max. load capacity	4 kg
Number of rollers	6
Roller size (length)	280 mm
Time setting range	-
Weight	4,5 kg



# Hermann



*Your professional manufacturing supplier  
of Laboratory Products*



Distributor in Argentina



**instrumentalia**

Gral. Urquiza 706 C1221ADF Bs. As. Argentina

Tel.: (54-11) 4932-3000 Fax: (54-11) 4931-6600

E-mail: [instrumentalia@pemint.com](mailto:instrumentalia@pemint.com)

<http://www.instrumentalia.com.ar>