

Hermann

Laboratory Products



*Your professional manufacturing supplier
of Laboratory Products*

- MAGNETIC HOTPLATE STIRRERS
- OVERHEAD STIRRERS
- VORTEX MIXERS
- ROTATORS
- TUBE ROLLER
- SHAKERS
- HOMOGENIZERS

» Shakers |

Hermann orbital and linear shakers are ideal for almost any vessel from microcentrifuge tubes through petri dishes and microtitre plates to conical flasks.

Los agitadores orbitales y lineales Hermann son ideales para casi cualquier recipiente de tubos de microcentrifuga a través de placas de Petri y placas de microtitulación a matraces cónicos.

HS-AJ-0660-PRO

LCD Digital Orbital Shaker



Performance Highlights

- Max. load capacity is 7.5kg with platform
- Double LCDs for independent display of time and speed
- Used with a wide variety of accessories
- Maintenance- free brushless DC motor
- Over speed detection and protection

Características Destacadas

- Max. capacidad de carga 7.5kg con plataforma
- Doble visor LCD para tiempo y velocidad
- Utilizado con una amplia variedad de accesorios
- Motor DC sin escobillas libre de mantenimiento
- Detección y protección de velocidad excesiva

HS-AJ-0360-W

Economical Digital Orbital Shaker



Performance Highlights

- Max. load capacity is 3kg with platform
- Speed range of 0~200 rpm
- A wide range of platforms

Características Destacadas

- Max. capacidad de carga 3kg con plataforma
- Rango de velocidad de 0 ~ 200 rpm
- Amplia gama de plataformas

Specifications	HS-AJ-0660-PRO	HS-AJ-0360-W
Voltage	100~240V,50/60Hz	100~240V,50/60Hz
Motor type	Brushless DC motor	DC motor
Orbital diameter	10 mm	20 mm
Timer display	LCD	-
Power	30W	20W
Speed range	100 - 500 rpm	0 - 200 rpm
Speed display	LCD	Scale
Time settings range	1 min - 19h 59 min	-
Weight	13.5 kg	3.1 kg
Dimension (WxDxH)	370×420×100mm	270×330×110mm
Permissible ambient temperature and humidity	5~40 °C, 80%	5~40 °C, 80%
Protection class	IP21	IP21
Max. load capacity (with platform)	7.5 kg	3 kg
Motor input power	28W	16W
Motor output power	15W	10W
Data connector	Rs232	-

Hermann



*Your professional manufacturing supplier
of Laboratory Products*



Distributor in Argentina



instrumentalia

Gral. Urquiza 706 C1221ADF Bs. As. Argentina

Tel.: (54-11) 4932-3000 Fax: (54-11) 4931-6600

E-mail: instrumentalia@pemint.com

<http://www.instrumentalia.com.ar>