

Psi precision system inc

5002 OSMETTE™ Automatic High Sensitivity Osmometer



Laboratory Model Freezing Point Osmometer for both 0.2 & 2.0 mL samples, featuring:

- Fully automatic operation over 0-3000 mOsm/kg range
- 30-60 sample/hour capability - Simple 2-step calibration
- 2 Ranges for higher accuracy
- Use with reusable or disposable tubes
- Probes last 2 to 4 times as long as comparable models

Comes with instruction manual, 24 2.0 mL sample tubes, 12 0.2 mL sample tubes, 100 & 500 mOsm/kg calibration standards, bath liquid, spare stirring wires, and a test tube rack. Power: 230v 50/60 Hz. OSMETTE A is shipped assembled and ready to use.